


FG GRUP

www.fggrup.com

REYA
HYBRID THROMBUS MANAGEMENT



REYA[®] THROMBECTOMY

www.biolas.net

The Reya Thrombectomy device is indicated for the mechanical opening of clotted vessels. With its tip design, it offers a unique treatment standard both mechanically and thrombolytically. The active tip gently disintegrates the thrombus. Thanks to its harmonic design, the active tip creates a vortex effect that does not allow the emboli to go distally.



The infusion port connected to the active tip allows TPA injection to the application area. Thus, the thrombus is easily dissolved by the mechano-chemical effect.

ADVANTAGES

- Preventing distal embolism
- Protecting the vascular wall
- Allowing TPA and contrast material infusion
- Excellent routing possibility
- Quick installation

REYA THROMBECTOMY SYSTEM	
REFERENCE NO	MODEL DIFFERENCES
ATA0007	90 cm Nitinol Shrink Wire
ATA0008	135 cm Nitinol Shrink Wire



Malıköy Başkent OSB Mah. Teknokent Blv. NO: 39
Sincan/ANKARA
Tel: +90 312 472 11 18/25 Fax: +90 312 472 11 26
info@fggrup.com.tr www.fggrup.com

www.fggrup.com